

HUAPLANT[®]

华普兰手术工具箱 HUAPLANT SURGICAL KIT

注册人/生产企业: 华磐(四川)医疗器械有限公司
住所/生产地址: 四川省资阳市雁江区外环路西三段222号7栋1楼1-6号
电话: 028-26933220
www.huaplant.com

Registered person/manufacturer:
Huaplant (Sichuan) Medical Equipment Co., Ltd
Address/Production Address:
1-6, 1st Floor, Building 7, No. 222, West Third Section, Outer Ring Road, Yanjiang District, Ziyang City, Sichuan Province
TEL: +86 28 26933220

Instruments.Surgical Kit

Basic Surgical KIT 基础手术工具箱

华普兰基础手术工具箱将允许您执行各种治疗程序。为确保最佳结果，请在使用工具包中包含的仪器和钻头时，遵循推荐的钻孔顺序。该套件包含了在标准的外科手术过程中所需的器械。


Huaplant basic surgical kit will allow you to perform various treatment procedures. To ensure optimal results, please follow the recommended drilling sequence when using the instruments and drill bits included in the toolkit. This kit includes the instruments required for standard surgical procedures.



基础手术工具箱 Basic Surgical KIT

Product	Type	Diameter	Cutting Length	Ref
	侧切钻 Lindemann Drill	Ø2.1	—	CLD21
	Ø2.25起始钻 Initial Drill	Ø2.25	10mm 14mm	CSID225 CSID225L
	扩孔钻 Reaming Drill	Ø2.8	10mm 14mm	PFD281 PFD281L
	扩孔钻 Reaming Drill	Ø3.4	10mm 14mm	CSFD3410 CSFD3414
	扩孔钻 Reaming Drill	Ø3.8	10mm 14mm	CSFD3810 CSFD3814
	扩孔钻 Reaming Drill	Ø4.1	10mm 14mm	CSFD4110 CSFD4114
	扩孔钻 Reaming Drill	Ø4.8	10mm 14mm	CSFD4810 CSFD4814
	皮质骨钻 Cortical Drill	Ø3.7	13.5mm	CSR37
	皮质骨钻 Cortical Drill	Ø4.1	13.5mm	CSR41
	皮质骨钻 Cortical Drill	Ø4.5	13.5mm	CSR45
	皮质骨钻 Cortical Drill	Ø5.2	13.5mm	CSR52

Product	Type	Diameter	Length	Ref
	植体螺丝起 Implant Driver	Hex 2.5	30mm	MDM08-25
	植体螺丝起 Implant Driver	Hex 2.5	35mm	MDM14-25
	植体螺丝起 Implant Driver	Hex 2.5	24mm	MDH08-25
	植体螺丝起 Implant Driver	Hex 2.5	30mm	MDH14-25
	六角螺丝起 Hex Driver	Hex 1.2	19mm	TDH13
	六角螺丝起 Hex Driver	Hex 1.2	25mm	TDH18
	平行杆 Parallel Pin	Ø2.2	18mm	PP2210
	平行杆 Parallel Pin	Ø2.2	24mm	PP2216
	延长杆 Drill Extension	Ø4.6	28.7mm	CSDE
	停止环 Stopper	Ø4.9	6mm	IRS90
	停止环 Stopper	Ø5.9	6mm	IRW90

Product	Type	Ref
	种植测量尺 Depth Gauge	HSDG02
	扭力扳手 Torque Wrench	TWI07

HF Submerged Surgical KIT

HF种植工具盒

华普兰HF种植手术工具盒，专为追求效率和节省空间设计的牙科医生提供的值得信赖的解决方案。该套件具有高质量的器械，可用于内六角连接种植体系统。根据指定的钻孔顺序方便地安排钻孔，优先考虑用户舒适度。

Huaplant HF Submerged Surgical Kit, a trustworthy solution designed specifically for dentists who pursue efficiency and space saving. This kit has high-quality instruments that can be used for internal hexagonal connection implant systems. Conveniently arrange drilling according to the specified drilling sequence, prioritizing user comfort.





种植工具盒 Submerged Surgical KIT

Product	Type	Diameter	Cutting Length	Ref
	引导钻 Guide Drill	Ø1.9	16.9mm	OPD2016
	球钻 Round Drill	Ø1.9	—	CSB18
	侧切钻 Lindemann Drill	Ø2.1	15.6mm	CLD21
	Ø2.2麻花钻 Twist Drill	Ø2.2	7mm	OID2207
	Ø2.2麻花钻 Twist Drill	Ø2.2	8.5mm	OID2208
	Ø2.2麻花钻 Twist Drill	Ø2.2	10mm	OID2210
	Ø2.2麻花钻 Twist Drill	Ø2.2	11.5mm	OID2211
	Ø2.2麻花钻 Twist Drill	Ø2.2	13mm	OID2213
	Ø3.0麻花钻 Twist Drill	Ø3.0	11.5mm	OID3011
	Ø3.0麻花钻 Twist Drill	Ø3.0	15mm	OID3015

Product	Type	Diameter	Cutting Length	Ref
	锥形钻 Taper Drill	F3.8	7mm	ORD3807
	锥形钻 Taper Drill	F3.8	8.5mm	ORD3808
	锥形钻 Taper Drill	F3.8	10mm	ORD3810
	锥形钻 Taper Drill	F3.8	11.5mm	ORD3811
	锥形钻 Taper Drill	F3.8	13mm	ORD3813
	锥形钻 Taper Drill	F4.2	7mm	ORD4207
	锥形钻 Taper Drill	F4.2	8.5mm	ORD4208
	锥形钻 Taper Drill	F4.2	10mm	ORD4210
	锥形钻 Taper Drill	F4.2	11.5mm	ORD4211
	锥形钻 Taper Drill	F4.2	13mm	ORD4213
	锥形钻 Taper Drill	F4.6	7mm	ORD4607
	锥形钻 Taper Drill	F4.6	8.5mm	ORD4608
	锥形钻 Taper Drill	F4.6	10mm	ORD4610
	锥形钻 Taper Drill	F4.6	11.5mm	ORD4611
	锥形钻 Taper Drill	F4.6	13mm	ORD4613
	锥形钻 Taper Drill	F5.0	7mm	ORD5007
	锥形钻 Taper Drill	F5.0	8.5mm	ORD5008
	锥形钻 Taper Drill	F5.0	10mm	ORD5010
	锥形钻 Taper Drill	F5.0	11.5mm	ORD5011
	锥形钻 Taper Drill	F5.0	13mm	ORD5013

Product	Type	Diameter	Length	Ref
	皮质骨钻 Cortical Drill	F3.8	36.2mm	OBD3808
	皮质骨钻 Cortical Drill	F4.2	36.2mm	OBD4208
	皮质骨钻 Cortical Drill	F4.6	36.2mm	OBD4608
	皮质骨钻 Cortical Drill	F5.0	36.2mm	OBD5008
	锥形平行杆 Taper Parallel Pin	Ø3.5	19mm	OPP3510
	锥形平行杆 Taper Parallel Pin	Ø4.0	19mm	OPP4010
	六角螺丝起 Hex Driver	Hex 1.2	23mm	CHD23
	六角螺丝起 Hex Driver	Hex 1.2	28mm	CHD28
	六角螺丝起 Hex Driver	Hex 1.2	19mm	SHHD19
	六角螺丝起 Hex Driver	Hex 1.2	26mm	SHHD26
	植体螺丝起 Implant Driver	Hex 2.1	30.7mm	FMD14NH
	植体螺丝起 Implant Driver	Hex 2.5	35.7mm	FMD14
	植体螺丝起 Implant Driver	Hex 2.1	25.5mm	SFHD2120
	植体螺丝起 Implant Driver	Hex 2.5	26.5mm	SFHD2520
	携带器螺丝起 Mount Driver	Ø4.8	19.9mm	CWM10
	携带器螺丝起 Mount Driver	Ø4.8	26.4mm	CWM15
	携带器螺丝起 Mount Driver	Ø4.8	15mm	SFHD4015
	携带器螺丝起 Mount Driver	Ø4.8	20mm	SFHD4020
	四方延长器 Torque Extension	S4.0	16mm	SHD16
	延长杆 Drill Extension	Ø4.5	27mm	DE0127

Product	Type	Ref
	种植测量尺 Depth Gauge	HSDG01
	扭力扳手 Torque Wrench	TWI07A

HX Submerged Surgical KIT

HX种植工具箱

华普兰HX种植手术工具箱，该套件采用便捷直观的设计方案，配备全套必要钻头工具，特别设计的超锋利钻头提供两种长度选择，通过最小化钻孔深度即可轻松穿透皮质骨，既优化了手术流程又显著提升了工作效率。

Huaplan HX Submerged Surgical Kit features a user-friendly and intuitive design, equipped with a complete set of essential drill tools. The specially designed ultra-sharp drills come in two length options, enabling effortless penetration through cortical bone by minimizing drilling depth. This optimizes surgical procedures while significantly enhancing work efficiency.



种植工具箱 Submerged Surgical KIT

Product	Type	Diameter	Cutting Length	Ref
	引导钻 Lance Drill	Ø2.3	15mm	HPD2315
	球钻 Round Drill	Ø1.9	—	CSB18
	侧切钻 Lindemann Drill	Ø2.1	15.6mm	CLD21
	Ø2.2起始钻 Initial Drill	Ø2.2	12mm	HID2210
	Ø2.2起始钻 Initial Drill	Ø2.2	15mm	HID2213
	扩孔钻 Reaming Drill	Ø2.8	12mm	HRD2810
	扩孔钻 Reaming Drill	Ø2.8	15mm	HRD2813
	扩孔钻 Reaming Drill	Ø3.2	12mm	HRD3210
	扩孔钻 Reaming Drill	Ø3.2	15mm	HRD3213

Product	Type	Diameter	Cutting Length	Ref
	扩孔钻 Reaming Drill	Ø3.6	12mm	HRD3610
	扩孔钻 Reaming Drill	Ø3.6	15mm	HRD3613
	扩孔钻 Reaming Drill	Ø4.0	12mm	HRD4010
	扩孔钻 Reaming Drill	Ø4.0	15mm	HRD4013
	扩孔钻 Reaming Drill	Ø4.4	12mm	HRD4410
	扩孔钻 Reaming Drill	Ø4.4	15mm	HRD4413
	扩孔钻 Reaming Drill	Ø4.8	12mm	HRD4810
	扩孔钻 Reaming Drill	Ø4.8	15mm	HRD4813
	扩孔钻 Reaming Drill	Ø5.2	12mm	HRD5210
	扩孔钻 Reaming Drill	Ø5.6	12mm	HRD5610

Product	Type	Diameter	Length	Ref
	停止环 Stopper	Ø4.9	6mm	IRS90
	停止环 Stopper	Ø4.9	4.5mm	IRS105
	停止环 Stopper	Ø5.9	6mm	IRW90
	停止环 Stopper	Ø5.9	4.5mm	IRW105
	平行杆 Parallel Pin	Ø2.2	18mm	PP2210
	平行杆 Parallel Pin	Ø2.2	24mm	PP2216
	延长杆 Drill Extension	Ø4.5	27mm	DE0127

Product	Type	Diameter	Length	Ref
	植体螺丝起 Fixture Driver	Hex 2.1	30.7mm	FMD14N
	植体螺丝起 Fixture Driver	Hex 2.1	30.7mm	FMD14NH
	植体螺丝起 Fixture Driver	Hex 2.5	29.7mm	FMD08
	植体螺丝起 Fixture Driver	Hex 2.5	35mm	FMD14
	植体螺丝起 Fixture Driver	Hex 2.5	24mm	FHD14
	植体螺丝起 Fixture Driver	Hex 2.5	30mm	FHD18
	六角螺丝起 Hex Driver	Hex 1.2	19.4mm	THD13
	六角螺丝起 Hex Driver	Hex 1.2	26.4mm	THD18
	携带器螺丝起 Mount Driver	Ø4.8	19.9mm	CWM10
	携带器螺丝起 Mount Driver	Ø4.8	26.4mm	CWM15

Product	Type	Ref
	扭力扳手 Torque Wrench	TWI07

Master KIT

大师工具箱

华普兰大师工具箱是集上颌窦内提、上颌窦外提、骨挤压、骨劈开、取骨、取植体、修整骨平面等多功能为一体的万能工具箱，安全实用操作简便。

Huaplant Master KIT is a versatile tool box that integrates functions such as maxillary sinus internal lifting, maxillary sinus external lifting, bone compression, bone splitting, bone extraction, implant removal, and bone plane correction. It is safe, practical, and easy to operate.



大师工具箱 Master KIT

Product	Type	Diameter	Ref
	球钻 Round Drill	Ø1.9	CSB18
	取骨钻 Trephine Drill	Ø3.0	TPE30
	内提钻 Internal Sinus Lift Drill	Ø2.9	CST28
	内提钻 Internal Sinus Lift Drill	Ø3.4	CST34
	内提钻 Internal Sinus Lift Drill	Ø4.0	CST40
	注水工具 Hydraulic Lifter	Ø4.5	THI45
	外开窗环钻 External Sinus Lift Drill	Ø6.0	CAC06
	外开窗球钻 External Sinus Lift Drill	Ø6.0	CAD06
	外开窗扩大器 External Sinus Lift Drill	Ø4.0	CAR06
	自体骨收集器 Bone Harvesting Drill	Ø4.5	TPP45
	菠萝钻 Bone Trimmer	Ø4.9	CPB05

Product	Type	Diameter	Ref
	骨锯 Round Bone saw	Ø7.0	CSS07
	骨锯 Round Bone saw	Ø9.0	CSS09
	骨挤压器 Bone Expander	Ø3.1	CSTS31
	骨挤压器 Bone Expander	Ø3.7	CSTS37
	骨挤压器 Bone Expander	Ø4.4	CSTS44
	骨挤压螺丝刀 Bone Expander Driver	—	MHDC2520
	骨挤压螺丝刀 Bone Expander Driver	—	MHDR2513


Product	Type	Diameter	Length	Ref
	停止环 Stopper	Ø4.9	12.0mm	IRS30
	停止环 Stopper	Ø4.9	11.0mm	IRS40
	停止环 Stopper	Ø4.9	10.0mm	IRS50
	停止环 Stopper	Ø4.9	9.0mm	IRS60
	停止环 Stopper	Ø4.9	8.0mm	IRS70
	停止环 Stopper	Ø4.9	7.0mm	IRS80
	停止环 Stopper	Ø7.0	7.5mm	ILS05
	停止环 Stopper	Ø7.0	7.0mm	ILS10
	停止环 Stopper	Ø7.0	6.5mm	ILS15
	停止环 Stopper	Ø7.0	6.0mm	ILS20
	停止环 Stopper	Ø7.0	5.5mm	ILS25
	停止环 Stopper	Ø7.0	5.0mm	ILS30

Product	Type	Ref
	粘膜剥离器 Sinus Elevator	LE01
	粘膜剥离器 Sinus Elevator	LE02
	粘膜剥离器 Sinus Elevator	LE03
	测量尺 Depth Gauge	LE04
	注水管 Irrigation Tube	GJG20
	扭力扳手 Torque Wrench	TWI07

GBR KIT
骨增量工具箱

垂直骨增量、穿翼、种植高阶 GBR 必备套装。梯度设计的骨凿与膜钉系统，让骨增量手术更精准高效，为种植体打造稳固“地基”。

A must-have Kit for vertical bone increment, wing penetration, and advanced GBR implantation. The gradient designed bone chisel and membrane nail system make bone augmentation surgery more precise and efficient, creating a stable "foundation" for the implant.





骨增量工具箱 GBR KIT

Product	Type	Diameter	Length	Ref
	骨钉 Bone Screw	Ø1.2	2.9mm	BS1203
			3.9mm	BS1204
			5.9mm	BS1206
			7.9mm	BS1208
	骨钉 Bone Screw	Ø1.4	3.9mm	BS1404
			5.9mm	BS1406
			7.9mm	BS1408
			9.9mm	BS1410
	骨钉 Bone Screw	Ø1.6	3.9mm	BS1604
			5.9mm	BS1606
			7.9mm	BS1608
			9.9mm	BS1610
			11.9mm	BS1612
	骨钉 Bone Screw	Ø1.8	3.9mm	BS1804
			5.9mm	BS1806
			7.9mm	BS1808
			9.9mm	BS1810
			11.9mm	BS1812

Product	Type	Diameter	Length	Ref
	帐篷钉 Tenting Screw	Ø1.6	6.3mm	TS1607
			7.3mm	TS1608
			8.3mm	TS1609
			9.3mm	TS1610
	膜钉 Membrane Tack	Ø0.9	3.45mm	MS0935
	膜钉 Membrane Tack	Ø0.9	3.95mm	MS0940

Product	Type	Diameter	Gingival Height	Ref
	基柱螺丝 Titanium Mesh Abutment	Ø2.0	1.0mm	LSC10
	基柱螺丝 Titanium Mesh Abutment	Ø2.0	2.0mm	LSC20
	基柱螺丝 Titanium Mesh Abutment	Ø2.0	3.0mm	LSC30
	基底帽 Healing Cap	Ø4.6	3.5mm	SC45

Product	Type	Diameter	Ref
	穿翼钻 Pterygoid Drill	Ø2.2	HIDW2230
	穿翼钻 Pterygoid Drill	Ø2.8	HIDW2830
	骨钉钻 Bone Screw Drill	Ø1.0	HIDN1010
	骨钉钻 Bone Screw Drill	Ø1.2	HIDN1210
	骨钉钻 Bone Screw Drill	Ø1.4	HIDN1410
	骨钉钻 Bone Screw Drill	Ø1.6	HIDN1610
	机用螺丝起 Machine Bone Screw Driver	Ø2.2	BMD21
	机用螺丝起 Machine Bone Screw Driver	Ø2.2	BMD31
	手用螺丝起 Handled Bone Screw Driver	Ø2.2	BHD64
	六角螺丝起 Hex Driver	Hex 1.2	THD18

Product	Type	Ref
	钛网 Titanium Mesh	HTMS20
	钛网 Titanium Mesh	HTMS30
	钛网 Titanium Mesh	HTMT20
	钛网 Titanium Mesh	HTMI30
	钛网 Titanium Mesh	HTMT30
	智能钛网(小) Titanium Mesh	HTMR15
	智能钛网(中) Titanium Mesh	HTMR16
	智能钛网(大) Titanium Mesh	HTMR18
	螺丝刀手柄 Screw Driver Handle	HD28
	取钉器 Tack Picker	TOMS
	持钉器 Tack Holder	CNH20
	牙骨锤 Dental Bone Hammer	BH160

SOS KIT
牙科维护工具盒

真正解决种植牙问题的万能工具盒。三大板块专攻种植机械并发症：

- 取断螺丝：可取市面上任意品牌的断螺丝（标配韩系，接受医生定制任意品牌接口）。
- 取种植体：两大方法任取多品牌失败、污染、断裂种植体（接受医生定制非标直径螺丝杆）。
- 取基台：两大方法任取断裂、损坏的空心、实心基台，以及滑扣的封闭螺丝、中央螺丝等。



The universal tool box that truly solves the problem of dental implants. Three major sectors specialize in complications of planting machinery:

- Remove broken screws: Any brand of broken screws on the market can be taken (standard Korean, accepting customized interfaces from doctors of any brand).
- Remove Implant: Two methods to retrieve failed, contaminated, or fractured implants from multiple brands (accepting customized non-standard diameter screw rods from doctors).
- Remove abutment: There are two methods to take the broken or damaged hollow or solid base, as well as the closed screw and central screw of the slider.

牙科维护工具盒 SOS KIT

Product	Type	Diameter	Ref
	钨钢钻 Reverse Guide Drill	Ø0.9	RSD09
	钨钢钻 Reverse Guide Drill	Ø1.1	RSD11
	螺丝取出器 Screw Remover	Ø0.9	RST09
	螺丝取出器 Screw Remover	Ø1.1	RST11
	取骨器 Trephine Drill	Ø4.8	CTD50R
	取骨器 Trephine Drill	Ø5.8	CTD60R
	取骨器 Trephine Drill	Ø6.8	CTD70R

Product	Type	Diameter	Ref
	限位套筒 Drill Guide	Hex 2.1	GD21
	限位套筒 Drill Guide	Hex 2.5	GD25
	植体连接器 Fixture Removal Screw	Ø1.6	FRS16
	植体连接器 Fixture Removal Screw	Ø2.0	FRS20
	植体取出器 Fixture Remover	Ø4.0	FRSD40
	植体取出器 Fixture Remover	Ø5.0	FRSD40
	连接器螺丝刀 Fixture Removal Screw Driver	—	HDH10
	一字螺丝刀 Slot Driver	—	AHD18
	基台取出器 Abutment Remover	—	TSHD23

Product	Type	Ref
	转接头 Adapter	MDT037
	固位扳手 Guide Holder	CW502
	扭力扳手 Torque Wrench	TWI07

Implant Cleaning Brush

种植体清洁刷

华普兰种植体清洁刷采用进口纯钛制造而成，与牙科种植手机连接后清洁种植体螺纹，对种植体周围的结石、菌斑、异物等进行清扫，可有效处理种植体周围炎，维护种植体牙周健康。

Huaplan Implant Cleaning Brush is made of imported pure titanium and is connected to a dental implant phone to clean the threads of the implant. It effectively removes stones, plaque, foreign objects, and other debris around the implant, effectively treating peri implantitis and maintaining periodontal health.



种植体清洁刷 Implant Cleaning Brush

Product	Type	Diameter	Length	Ref
	常规钛刷 Titanium Brush	Ø1.2	30.6mm	ICB3012
	常规钛刷 Titanium Brush	Ø1.5	30.6mm	ICB3015
	常规钛刷 Titanium Brush	Ø1.2	36.1mm	ICB3612
	常规钛刷 Titanium Brush	Ø1.5	36.1mm	ICB3615
	T型钛刷 T-shaped Titanium Brush	Ø1.8	32mm	TCB3218

Bone Flattening KIT

骨修整工具盒

华普兰骨修整工具盒为种植手术中出现基台就位不良的情况提供了精准解决方案。骨修整工具盒可与多款品牌的种植体系统配套使用，不仅能有效防止骨屑飞溅伤到植体内部，还能预留有效骨量，更最大限度保留牙槽嵴。

Huaplan bone flattening kit provides a precise solution for cases of poor abutment positioning during implant surgery. Bone flattening kit can be used in conjunction with multiple brands of implant systems, which not only effectively prevents bone debris from splashing into the interior of the implant, but also reserves effective bone mass and maximizes the preservation of the alveolar ridge.



骨修整工具盒 BLG KIT

Product	Type	Ref
	引导柱(BLT NC) Bone Flattening Guide Pin for Straumann® BLT NC	BLG01
	引导柱(BLT RC) Bone Flattening Guide Pin for Straumann® BLT RC	BLG02
	引导柱(诺保科 NP) Bone Flattening Guide Pin for Nobel Biocare™ NP	BLG03
	引导柱(诺保科 RP) Bone Flattening Guide Pin for Nobel Biocare™ RP	BLG04
	引导柱(奥齿泰 M) Bone Flattening Guide Pin for Osstem® M	BLG05
	引导柱(奥齿泰 R) Bone Flattening Guide Pin for Osstem® R	BLG06
	引导柱(ADIN OS) Bone Flattening Guide Pin for ADIN® OS	BLG07
	引导柱(登腾) Bone Flattening Guide Pin for Dentium®	BLG08

* 6 bone flattening guide pins can be freely selected according to needs

Product	Type	Diameter	Ref
	平骨钻 Bone Flattening Drill	Ø4.2	BLD42
	平骨钻 Bone Flattening Drill	Ø4.8	BLD48
	平骨钻 Bone Flattening Drill	Ø5.4	BLD54
	平骨钻 Bone Flattening Drill	Ø6.0	BLD60
	六角螺丝起 Hex Driver	Hex 1.27	HHD22

Zygomatic Implant Tool

穿颧种植工具

- 球钻：穿颧种植手术中用于定位种植体的位置并确定间距
- 穿颧钻：用于穿颧种植备洞
- 穿颧测量尺：探诊种植备洞的窝洞深度
- Round Drill: Used in zygomatic implant surgery to locate the position of the implant and determine the spacing
- Zygomatic Drill: Used for preparing holes for zygomatic implantation
- Zygomatic Depth Gauge: Probing the depth of the osteotomy during implant site preparation

穿颧种植工具 Zygomatic Implant Tool

Product	Type	Diameter	Length	Ref
	长球钻 Long Round Drill	Ø2.9	70mm	CSBL29
	穿颧钻 Zygomatic Drill	Ø2.9	95mm	HIDW2840
	穿颧钻 Zygomatic Drill	Ø3.5	100mm	HIDW3540
	穿颧测量尺 Zygomatic Depth Gauge	Ø2.1	167mm	HSDG04

Pterygoid Implant Tool

穿翼种植工具

- 植入扳手：用于拧紧或旋出种植体
- 穿翼锥：在植入位点挤压扩张备孔通道
- 穿翼钻：用于穿翼种植备洞
- Handled Fixture Driver: Used for tightening or unscrewing implants
- Pterygoid Awl: Squeezing and expanding the pore preparation channel at the planting site
- Pterygoid Drill: Used for preparing holes for pterygoid implant surgery

穿翼种植工具 Pterygoid Implant Tools

Product	Type	Diameter	Length	Ref
	植入扳手 Handled Fixture Driver	Hex 2.5	138mm	TPP1
	植入扳手 Handled Fixture Driver	Hex 2.4	138mm	TPP3
	穿翼锥 Pterygoid Awl	Ø2.1	162mm	TPP22
	穿翼锥 Pterygoid Awl	Ø2.7	162mm	TPP28
	穿翼钻 Pterygoid Drill	Ø2.2	45mm	HIDW2230
	穿翼钻 Pterygoid Drill	Ø2.8	45mm	HIDW2830